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Shenzhen Lvshiyuan Biotechnology Co., Ltd, was founded in 2003 by a senior scholar who had studied in America for many years. It is a high-tech enterprise specializing in the development and production of the most advanced In vitro diagnostic reagents(IVD), food safety test reagent, animal disease diagnostic reagent, small molecule antigen & antibody, committed to the industrialization of biological reagent around the world.

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LSY get the ISO9001 certificate of quality management system, The review of ISO13485 has been passed, we can get this certificate soon. we has achieved many qualification certificates, such as National HI-TECH enterprise certificate, Shenzhen HI-TECH enterprise certificate, Postdoctoral innovation practice base, AAA-grade credit enterprise, Medical Equipment Production License etc. LSY also get a number of patent certificate.

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In the 21st biological high-tech era, LSY takes the “Health care, benefit to society; Focus on products, lead the industry; Create benefits, repay shareholders; perfect mechanism, succeed employees. “as our Enterprise mission and adheres to the Quality philosophy of “Quality makes brand, brand is life.”, providing the advanced biological technology products and the perfect technical service for the customers.

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Guardian Independent Certification Ltd

Registered in England

Sovereign House, 212-224 Shaftesbury Avenue London England WC2H 8HQ

Accredited by Member of IAF MLA

IAF-ANZ registration no. 53510508UK, www.jas-anz.org/register

ISO 13485 认证证书 ISO 13485 certification

Certificate CN20/42084

The management system of

Shenzhen Lvshiyuan Biotechnology Co., Ltd.

101, 201, 301, Building D, No. 2, Industrial Avenue, Buxin Village,
Buxin Community, Dapeng Subdistrict Office, Dapeng New District,
Shenzhen, Guangdong, 518120, P.R. China

has been assessed and certified as meeting the requirements of

ISO 13485:2016 EN ISO 13485:2016

For the following activities

**Design, Manufacture and Distribution of Dry Fluorescent
Immunoassay Instrument and In Vitro Diagnostic Test Kits
(ELISA, Colloidal Gold) for SARS-CoV-2, Influenza A and Influenza B.**

This certificate is valid from 2 March 2021 until 13 June 2023
and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 26 May 2023
Issue 2. Certified since 14 June 2020

Authorised by

SGS United Kingdom Ltd
Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK
t +44 (0)151 350-6666 f +44 (0)151 350-6600 www.sgs.com

HC SGS 13485 2016 0118

Page 1 of 1

This document is issued by the Company subject to its General Conditions of Certification Services accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <http://www.sgs.com/en/certified-clients-and-products/certified-client-directory>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

对外贸易经营者备案登记表 Foreign trade operator registration form

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报关注册登记证书 Customs registration certificate

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医疗器械生产许可证 Medical device production license

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Obtained screenshots of foreign standard certification

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1 von 1

12.02.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with $CT < 25$, 23 samples with CT between 25 and 30, and 9 samples with $CT > 30$. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 µL of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that many other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

Contact

Email: sarscov2ivd@pei

Last updated: 12.02.2021

Overview of SARS-CoV-2 Antigen Rapid Tests Assessed as Reflecting the Current State of the Art

Name of Test	Manufacturer (Distributor)
Panbio™ COVID-19 Ag Rapid Test Device (NASOPHARYNGEAL)	Abbott Rapid Diagnostics Jena GmbH
RIDA®QUICK SARS-CoV-2 Antigen	R-Biopharm AG
SARS-CoV-2 Rapid Antigen Test	SD BIOSENSOR (Roche Diagnostics GmbH)
NADAL® COVID-19 Ag Schnelltest	nal von minden gmbh
STANDARD™ F COVID-19 Ag FIA	SD BIOSENSOR
STANDARD™ Q COVID-19 Ag Test	SD BIOSENSOR
BIOSYNEX COVID-19 Ag BSS	BIOSYNEX SWISS SA
MEDsan® SARS-CoV-2 Antigen Rapid Test	MEDsan GmbH
TestNOW® - COVID-19 Antigen	Affimedix
NowCheck® COVID-19 Ag Test	BIONOTE
Coronavirus Ag Rapid Test Cassette (Swab)	Zhejiang Orient Gene Biotech Co.,Ltd
Sofia SARS Antigen FIA	Quidel Corporation
COVID-19 Ag Test Kit	Guangdong Wesail Biotech Co., Ltd.
CLINITEST® Rapid COVID-19 Antigen Test	Siemens Healthineers
ESPLINE® SARS-CoV-2	Fujirebio Inc. (Mast Diagnostica GmbH)
BD Veritor™ System for Rapid Detection of SARS-CoV-2	Becton Dickinson
GenBody COVID-19 Ag	IVC Pragen Healthcare
LumiraDx SARS-CoV-2 Ag Test	LumiraDX
Exdia COVID-19-Ag-Test	Precision Biosensor Inc. (Axon Lab AG)
SARS-CoV-2 Ag Rapid Test (FIA)	Wantai (Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.)
SARS-CoV-2 Antigen Schnelltest	Xiamen Boson Biotech Co., Ltd (Medicovid-AG; technomed GmbH; Löwe Medizintechnik)
COVID-19 Antigen Schnelltest (Colloidal Gold)	Joinstar Biomedical Technology Co., Ltd (CIV care impuls Vertrieb)
mö-screen Corona Antigen Test	Mölab GmbH
Rapid SARS-CoV-2 Antigen Test Card	MP Biomedicals Germany GmbH
Lyher Novel Coronavirus (COVID-19) Antigen Test Kit (Colloidal Gold)	Hangzhou Laihe Biotech Co., Ltd. (Lissner Qi GmbH)
AMP Rapid Test SARS-CoV-2 Ag	Ameda Labordiagnostik GmbH
Clungene COVID-19 Antigen Rapid Test	Hangzhou Clongene Biotech Co., Ltd.
Gensure™ COVID-19 Antigen Rapid Test Kit	GenSure Biotech Inc.
SARS-CoV-2 Antigen Rapid Test Kit	Beijing Lepu Medical Technology Co., Ltd
Hightop SARS-CoV-2 (Covid-19) Antigen Rapid Test	Qingdao Hightop Biotech Co., Ltd.
Rapid Covid-19 Antigen Test (Colloidal Gold)	Anbio (Xiamen) Biotechnology Co., Ltd

Name of Test	Manufacturer (Distributor)
Safecare COVID-19 Ag Rapid Test Kit (Swab)	Safecare Biotech Hangzhou Co., Ltd.
QuickProfile Covid-19 Antigen Test Card	LumiQuick Diagnostics, Inc.
Covid 19 Antigen Schnelltest	BioRepair GmbH
Green Spring SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	Shenzhen Lvshiyuan Biotechnology Co., Ltd.
CAT Antigen Covid Rapid Test	Oncosem Onkolojik Sistemler San. Ve Tic. A.S.
ScheBo SARS-CoV-2 Quick Antigen	ScheBo Biotech AG
Nova Test SARS-CoV-2 Antigen Rapid Test Kit	Atlas Link Technology Co., Ltd.
Toda Coronadiag Ag	Toda Pharma
Humasis COVID-19 Ag Test	Humasis Co., Ltd.
Beijing Hotgen Biotech Co., Ltd.	Neuartiges Coronavirus (2019-nCoV)- Antigentest (Kolloidales Gold); Novel Coronavirus 2019-nCoV Antigen Test (Colloidal gold)
Xiamen AmonMed Biotechnology Co.,Ltd.	COVID-19 Antigen Rapid Test Kit (Colloidal Gold)
Canea COVID-19 Antigen Schnelltest	Core Technology Co., Ltd.
fluorecare COVID-19 SARS-CoV-2 Spike Protein Test Kit (Colloidal Gold Chromatographic Immunoassay)	Shenzhen Microprofit Biotech Co., Ltd
Tetsealabs® Rapid Test Kit COVID-19 Antigen Test Cassette	Hangzhou Testsea Biotechnology Co., Ltd
Lysun COVID-19 Antigen Rapid Test Device (Colloidal Gold)	Hangzhou Lysun Biotechnology Co., Ltd.

欧盟通用和互认清单 On the EU common list & Mutual Recognition list

https://covid-19-diagnostics.jrc.ec.europa.eu/devices/detail/2109

COVID-19 In Vitro Diagnostic Devices and Test Methods Database

Home > COVID-19 In Vitro Diagnostic Medical Devices > COVID-19 In Vitro Diagnostic Medical Device - detail

COVID-19 In Vitro Diagnostic Medical Device - detail

< Previous Next >

Green Spring SARS-CoV-2 Antigen-Rapid test-Set

Manufactured by Shenzhen Lvshiyuan Biotechnology Co., Ltd., China - <https://www.lsybt.com/>

Device identification number	2109
CE Marking	Yes
HSC common list	Yes
HSC mutual recognition	Yes
Format	Manual, Near POC / POC
Physical Support	Lateral flow
Target	Antigen
Specimen	Anterior nasal swab, Nasal swab, Nasopharyngeal swab, Oropharyngeal swab
Pathogens detected	Coronaviruses (HCoV), SARS-CoV
Lineages detected	B.1.1.7 (United Kingdom) , B.1.351 (South Africa) , P.1 (Japan/Brazil)
Commercial Status	Commercialised
Last Update	2021-07-07 05:24:52 CET
Comments	Der Green Spring® SARS-CoV-2-Antigen-Schnelltest dient dem schnellen qualitativen Nachweis des Nukleocapsid-Protein-Antigens von SARS-CoV-2 in menschlichen Nasen-, Nasen-Rachen oder Rachenabstrichproben. Die Ergebnisse dienen dem Nachweis von SARS-CoV-2-Antigenen.

比利时白名单 Whitelist of Belgium

Manufacturer	Device Name	CE Marking	HSC Common List	HSC Mutual Recognition	Format	Physical Support	Target	Specimen	Pathogens detected	Lineages detected	Commercial Status	Last Update	Comments	Whitelist Status
Liclear biotech (Hangzhou)	SARS-CoV-2 Nucleocapsid (N) Antigen Rapid Test Cassette	95.1	99.3				NP swab							No
LumiQuick Diagnostics	QuickProfile COVID-19 Antigen Test Strip	94.0	99.0				NP swab							Yes
MEDsan	SARS-CoV-2 Antigen Rapid Test	92.5	99.8				Nasal swab/OP swab							Yes
MP Biomedicals	Rapid SARS-CoV-2 Antigen Test Card	96.4	99.0				NP swab/OP swab							Yes
Multi-G	COVID19CHECK-GEN	92.5	99.2				NP swab/OP swab							No
	COVID19CHECK-NAS (COVID-19 antigen rapid test cassette)	98.1	99.2				Nasal swab							No
	COVID19CHECK-SAL (Sars-Cov-2 Antigen Rapid Test)	96.8 (Cts 33)	99.1				Saliva							No
MyLab Discovery Solutions Pvt	PathoCatch COVID-19 Antigen Lateral Flow Test Device	92.0	100.0				Nasal swab							No
nal von minden	NADAL COVID-19 Ag Rapid Test	97.6	99.9				NP swab/OP swab							Yes
Nanjing Vazyme Medical Technology	Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) Antigen Detection Kit (Colloidal Gold-Based)	97.6	99.3				Nasal swab/OP swab							No
NanoEntek	COVID-19 Antigen Saliva Test Kit (Colloidal Gold) (Cassette)	94.7	100.0				NP swab							Yes
	COVID-19 Antigen Test Kit (Colloidal Gold)	98.1	99.3				Saliva							No
Nantong Diagnosis Biotechnology	COVID-19 Antigen Detection Kit	98.8/98.5	100.0				NP swab/OP swab							No
New Gene (Hangzhou) Bioengineering	COVID-19 Antigen Detection Kit	97.3/95.7/95.1	99/99/99.1				Nasal swab/OP swab/Sputum							Yes
ONCOSEM Onkoloogik Sistimler San. ve Tic. A.Ş.	2019-nCoV Antigen Rapid Test Kit	97.6	99.3				NP swab							No
Ortho-Clinical Diagnostics	VITROS SARS-CoV-2 Antigen test	97.8	99.2				NP swab							Lab test
PCL	PCL COVID19 Ag Gold	90.8 (Ct<30)/91.7	99.5/100				Saliva/NP swab							No
Prestige Diagnostics	COVID-19 nCoV Antigen Device	90.9	99.1				NP swab							No
PRIMA Lab	PRIMACOVID COVID-19 Antigen Rapid Test	96.4/92.9	99.2/100				NP swab/Nasal swab							Yes
Qinadao Hightop Biotech	SARS-CoV-2 Antigen Rapid Test (Immunochromatography)	92.7/95.0/95.0	99.8/99.8/99.8				Nasal swab/NP swab/OP swab							Yes
Quidel	Sofia SARS Antigen FIA	96.7	100.0				NP swab/Nasal swab							Yes
	Sofia 2 Flu + SARS Antigen FIA	95.2	100.0				NP swab/Nasal swab							No
Roche Diagnostics	Elecsys SARS-CoV-2 Antigen	94.5	99.9				NP swab/OP swab							Lab test
SD Biosensor (distributed by Roche)	SARS-CoV-2 Rapid Antigen Test	96.5	99.7				NP swab							Yes
Shenzhen Huan Biosci Technology	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	94.6	99.1				Nasal swab							No
Shenzhen Landwind Biotechnology	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	95.7	99.2				Nasal swab							No
Shenzhen Lvshiyuan Biotechnology	Green Spring SARS-CoV-2 Antigen-Rapid Test Kit (Colloidal Gold)	96.8/98	100.0				Nasal swab/NP swab							Yes
Shenzhen Microprofit Biotech	Fluorecare SARS-CoV-2 Test Kit	98.8/100	100.0				NP swab/Nasal swab							No
	Fluorecare SARS-CoV-2 Spike Protein Test Kit	96.4	100.0				NP swab							No
	SARS-CoV-2 Antigen Test Kit (SC0201 + SC0202)	92/100/96	100.0				NP swab/OP swab/Saliva							No
Shenzhen Ultra Diagnostics Biotech	SARS-CoV-2 Antigen Test kit (SC0203 + SC0204)	92/95.7/97.3	100/99/99				NP swab/OP swab/Sputum							No
Shenzhen YHLO Biotech Co	GINE-2019-nCoV Ag	96.5/97.4	99.3/99.3				NP swab/Nasal swab							Yes
Surescreen	COVID-19 Coronavirus Rapid Antigen Test Cassette	93.3	100.0				NP swab/OP swab							No
Todapharma	TODA Coronadiag Ag	96.6	100.0				NP swab/OP swab							Yes
Ulti med Products	COVID-19 Antigen Test (Nasopharyngeal Swab)	96.4	99.2				NP swab							No
	COVID-19 and Influenza A+B Antigen Combo Rapid Test (Nasopharyngeal Swab)	96.4	99.2				NP swab							No
Van Oostveen Medical	Coronavirus Ag Rapid Test Cassette (Swab)	96.7	99.2				NP swab							No
VivaDiag Pro	VivaDiag Pro SARS-CoV-2 Ag Rapid Test	96.1/96.1/97.0	100.0				NP swab/OP swab/Nasal swab							Yes
VivaChek Biotech	VivaDiag SARS-CoV-2 Ag Rapid Test	95.0	100.0				NP swab/OP swab/Nasal swab							No
	VivaDiag SARS-CoV-2 Ag Saliva Rapid Test	98.3	100.0				Saliva							No

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Green Spring ®SARS-CoV-2 Antigen Rapid Test Kit

Inverkehrbringer

Shenzhen Lvshiyuan Biotechnology Co.,Ltd
D Building, National Biological Industrial Park of
Marinelife, No.2 Binhai Road, Dapeng, Shenzhen,
China

Hersteller

Shenzhen Lvshiyuan Biotechnology Co.,Ltd D
Building, National Biological Industrial Park of
Marinelife, No.2 Binhai Road, Dapeng, Shenzhen,
China

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英国注册证书 Registration Certificate of UK

Medicines & Healthcare products
Regulatory Agency

Medicines & Healthcare products
Regulatory Agency

10 South Colonnade
Canary Wharf
London
E14 4PU
United Kingdom

+44 (0) 20 3080 6000
gov.uk/mhra

PureUKCA Ltd
59
St. Martin's Lane
Middlesex
London
WC2N 4JS
England, United Kingdom

09 September 2021

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We are pleased to confirm that the application to register or update an existing registration for the following manufacturer, which you submitted on **09 September 2021** has been reviewed:

Application reference: **2021090901215399**

Manufacturer organisation: **Shenzhen Lvshiyuan Biotechnology Co.,Ltd.**

Address:

**D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng
New District, Shenzhen, 518120 China
Shenzhen
518120
China**

Manufacturer registration status: **Registered**

Device(s):

GMDN term	Status	Comment
SARS-CoV-2 antigen IVD, kit, immunochromatographic test (ICT), rapid	Registered	

Please note this letter **does not** represent any form of accreditation, certification or approval by the UK Competent Authority.

If you stop placing devices on the market or if you are not complying with the Regulations, you should inform us so that we can amend our records. You should be aware that it is an offence to place on the market UKCA or CE marked devices that do not comply with the regulations.

Please inform us of the following chargeable changes:

1. **company/organisation information e.g. name and address**
2. **additional devices (GMDN code or term)**

Please also use the Devices Online Registration Database (DORS) to tell us of the following changes e.g. removal/discontinuation of a device (GMDN) or product from your registration record, change of contact person, telephone number and/or email address, for which payment of our statutory fee does not apply.

Please note that the name and address of manufacturer, UK Responsible Person or Authorised Representative (Northern Ireland only) and devices that have been registered will be published on our [Public Access Registration Database \(PAR\)](#).

The account number for your company/organisation is **0000018481**.

Yours sincerely,

Ngozi Onyeukwu
 Device registrations service
 Devices division
 MHRA

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MINISTRY OF HEALTH

PHARMACY AND POISONS BOARD

(Section 3B(2)(e) of the Pharmacy and Poisons Act, Cap 244 Laws of Kenya)

IN-VITRO DIAGNOSTIC EMERGENCY USE AUTHORIZATION

This Emergency Use Authorization is issued to **Shenzhen Lvshiyuan Biotechnology Co., Ltd.** for distribution and sale of **SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)**

Emergency use Authorization (EUA) No.	MD/2021/10326
EUA valid until	End of COVID -19 Pandemic or EUA revocation
Device category	Medical Device class C/D
GMDN	N/A
GMDN Term	N/A
Intended purpose	Qualitative detection of SARS-CoV-2 nucleocapsid antigen (Ag) in human nasal swab specimens from individuals who are suspected of COVID-19 by use of immunochromatographic technique
Conditional Approval	For Screening Purposes Only

EUA No.: MD/2021/10326

Date of Authorization: September 15th 2021

MINISTRY OF HEALTH

THE PHARMACY AND POISONS ACT
(Cap. 244, Sub. Leg.)
(The Pharmacy and Poisons Rules)

WHOLESALE DEALER'S LICENCE

Messrs..... **BELLWAY PHARMACEUTICALS LIMITED -RIARA RD** of **JAIKULPHARM2020@YAHOO.COM**

carrying on business at **330/192** are hereby authorized to
sell poisons by way of wholesale dealing.

23-02-2021

Date

Note. (i) This licence expires on 31st day of December **2021**

Fee: KShs. 30000

Licence No. BU202108744

MINISTRY OF HEALTH
PHARMACY AND POISONS BOARD

Telegram: "MINHEALTH" Nairobi
Telephone: 020-2716905/6, 020-3562107
Cellphone: 0733-884411/0720 608811
Fax: 2713409
Email: admin@pharmacyboardkenya.org
Website: www.pharmacyboardkenya.org

Pharmacy & Poisons Board Head
Along Lenana Road
P. O. Box 27663-00506
NAIROBI

When replying please quote our ref No:

PPB/PER/GEN/VOL.I/21/35

Date: September, 15th 2021

MAH (MANUFACTURER/SUPPLIER)

Shenzhen Lvshiyuan Biotechnology Co., Ltd,
101,201,301 Building, No.2
Industrial Avenue, Buxin Village,
Buxin Community Dapeng Subdistrict Office,
Dapang New District Shenzhen,
Guangdong P.R China.

Through'

Local Technical Representative

Bellway Pharmaceutical Limited,
330/192 Riara Road, Kilimani,
P.O. Box 25443-00603,
Nairobi

Dear Sir/Madam,

**RE: APPLICATION FOR IN-VITRO DIAGNOSTIC MEDICAL DEVICE
EMERGENCY USE AUTHORIZATION OF SARS-CoV-2 Antigen Rapid Test Kit
(Colloidal Gold)**

Reference is made to the above subject.

The Pharmacy and Poisons Act (Cap 244) provides for the regulatory oversight of Medical Devices including In-Vitro Diagnostics (IVDs) in Kenya by the Pharmacy and Poisons Board "the Board". Please note that the Expert committee on medicals devices & IVDs of the Board is tasked with reviewing applications for Covid-19) test kits. The Board's policy requires that all IVDs meet certain conditions prior to approval; including submission of product dossier and validation protocols and reports that evidence threshold of key performance parameters. The validating laboratory must be an accredited laboratory, such as, Kenya Medical Research Institute that is recognized by the Board.

荷兰注册证 Registration Certificate in Netherlands

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Ministry of Health, Welfare and Sport
CIBG
P.O. Box 16114
2500 BC The Hague
THE NETHERLANDS

STATEMENT

The undersigned herewith declares that according to the Decree on In-Vitro Diagnostics, which is based on the European Directive 98/79/EC concerning in-vitro diagnostic medical devices,

Kingsmead Service B.V.
Zonnehof 36
2632 BE Nootdorp
THE NETHERLANDS

acts as authorised representative of the manufacturer.

The manufacturer:

Shenzhen Lvshiyuan Biotechnology Co.,Ltd
D Building, National Biological Industrial Park of Marinelife, No.2 Binhai Road, Dapeng
Shenzhen 518120
CHINA

is authorised to manufacture and/or supply the medical device/devices mentioned below:

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

This device/these devices may be placed on the Dutch market and on the markets of the other Member States of the European Union, and be exported to States and territories outside the European Union. This free sale certificate may only be used for export outside the European Union.

The present statement is drawn up at the request of the interested party in order to be submitted to the Health Authorities of SAUDI ARABIA .

This statement is valid until May 26, 2022.

The Hague, January 10, 2022

On behalf of State Secretary for Health, Welfare and Sport
Farmatec | CIBG

Dr. M.J. van de Velde
Mr. M.J. van de Velde
Head of Department

Our reference: 20217418
Certificate number: 34659

Ministry of Health, Welfare and Sport
CIBG
P.O. Box 16114
2500 BC The Hague
THE NETHERLANDS

STATEMENT

The undersigned herewith declares that according to the Decree on In-Vitro Diagnostics, which is based on the European Directive 98/79/EC concerning in-vitro diagnostic medical devices,

Kingsmead Service B.V.
Zonnehof 36
2632 BE Nootdorp
THE NETHERLANDS

acts as authorised representative of the manufacturer.

The manufacturer:

Shenzhen Lvshiyuan Biotechnology Co.,Ltd
D Building, National Biological Industrial Park of Marinelife, No.2 Binhai Road, Dapeng
Shenzhen 518120
CHINA

is authorised to manufacture and/or supply the medical device/devices mentioned below:

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

This device/these devices may be placed on the Dutch market and on the markets of the other Member States of the European Union, and be exported to States and territories outside the European Union. This free sale certificate may only be used for export outside the European Union.

The present statement is drawn up at the request of the interested party in order to be submitted to the Health Authorities of UNITED ARAB EMIRATES.

This statement is valid until May 26, 2022.

The Hague, January 10, 2022

On behalf of State Secretary for Health, Welfare and Sport
Farmatec | CIBG

Dr. M.J. van de Velde
Mr. M.J. van de Velde
Head of Department

Our reference: 20217418
Certificate number: 34660

马来西亚自测证 Self-test approval of MDA Malaysia

PIHAK BERKUASA PERANTI PERUBATAN
Medical Device Authority
KEMENTERIAN KESIHATAN MALAYSIA
Ministry of Health Malaysia
Aras 6, Prima 9, Prima Avenue II,
Blok 3547, Persiaran Apec,
63000 Cyberjaya, Selangor
Malaysia.

Ref : (76) MDA.600-1/6/27

Date : 19 November 2021

SPD Scientific (M) Sdn Bhd
F-G-45, Block F, Jalan PJU 1A/3
Taipan 2 Damansara, Ara Damansara
47301 Petaling Jaya, Selangor Darul Ehsan
(Attention to: Hor Zian Khang)

Dear Sir/Madam,

CONDITIONAL APPROVAL FOR IMPORTATION AND DISTRIBUTION OF MEDICAL DEVICE (COVID-19 SELF TEST KIT)

With reference to the above, I wish to inform that the Authority grants your establishment a conditional approval for the importation and distribution of medical device as listed in **Appendix 1**.

2. Please be informed that the validity of this conditional approval is from **19/11/2021** to **19/11/2022** and is subject to the following:

- i) Your establishment shall ensure that the medical device under this conditional approval complies with safety and performance requirements as stipulated in Medical Device Act 2012 (Act 737);
- ii) Your establishment shall adhere to the conditions as stipulated in **Appendix 2**.
- iii) The use of COVID-19 self test kit shall be limited for screening purpose only and all test results need further confirmation using RT-PCR.

3. This conditional approval for importation and distribution of this medical device is an interim measure in response to the current public health need during COVID-19 pandemic. This letter shall not be used for the purpose of promoting or advertising of the product and it does not exempt you from abiding to any laws or requirements by any other authorities of Malaysia.

Thank you,

(AHMAD SHARIFF BIN HAMBALI)
Chief Executive
Medical Device Authority
Ministry of Health Malaysia

Appendix 1

Medical Device Details


Name of Medical Device	: Green Spring SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)
Brand/Model	: Green Spring
Identifier	: GF102B1S
Sample type	: Nasal swab
Intended Use	: It is used to qualitatively detect SARS-CoV-2 nucleocapsid protein antigen in human nasal swab specimen. It is meant for self-testing in home
Brief Description	: The Green Spring® SARS-CoV-2 Antigen Rapid Test is a qualitative, membrane based immunoassay for the detection of SARS-CoV-2 nucleocapsid protein antigens. The test line region is coated with SARS-CoV-2 antibody. The sample reacts with the SARS-CoV-2 antibody in the test line region. If the specimen contains SARS-CoV-2 antigens, a colored line appears in the test line region as a relevant result. As a procedural control, a colored line appears in the control line area indicating that the correct volume of sample has been added and membrane wetting has occurred correctly. Result can be read in 15-20 minutes
Lot Number	: 20210401
Manufacturer's name	: Shenzhen Lvshiyuan Biotechnology Co., Ltd.P.R.C

Appendix 2

Conditions:

- 1) The conditional approval for importation and distribution of the medical device listed in Appendix 1 is valid for one year.
- 2) An establishment given the conditional approval shall—
 - i) collect data related to safety and performance of the medical device and shall submit the report to the Authority on a regular basis or when it is required by the Authority;
 - ii) submit any information requested by the Authority within the prescribed period;
 - iii) comply with any directions issued by the Authority from time to time;
 - iv) comply with labelling requirements stipulated in Sixth Schedule of the Medical Device Regulations 2012, in particular instruction for use and disposal method, including using infographic, to make it easily understood by lay persons;
 - v) provide suitable and adequate storage to ensure proper conservation of the medical device in accordance with the manufacturer's instruction;
 - vi) perform inspection on the primary packages of the medical device and any breached packages shall be disposed off appropriately;
 - vii) distribute the medical device only to licensed community pharmacies and healthcare institutions and they may sell the medical device online subject to appropriate distribution method specified by the manufacturer;
 - viii) establish adequate precautions and control to prevent deterioration or damage of the medical devices up until the point of use;
 - ix) ensure the delivery of medical devices adhere to the conditions specified by the manufacturer;
 - x) provide documentation of all medical devices supplied to customers, the quantity supplied, the batch or lot number and/or model and serial number;
 - xi) establish and maintain an appropriate distribution records up to retail distribution of the medical device to the end-user;
 - xii) keep the record of delivery transactions as the proof of supply of the medical device to customers;
 - xiii) dispose off medical device in accordance with regulatory requirements and any other applicable statutory requirements; and
 - xiv) not carry out any secondary assembly activities on the medical device unless the manufacturers instruction states otherwise;
- 3) All information pertaining to this medical device including all supporting documents shall be kept at the premises and shall be made available upon request by the Authority.
- 4) An establishment shall establish and maintain a post-market surveillance system to monitor the traceability of the medical device throughout the supply chain.
- 5) The Authority reserves the right to make a visit or inspection to the establishment at any time without prior notice.
- 6) The Authority may revoke the conditional approval or may take legal action should the establishment fails to comply with any conditions imposed by the Authority.
- 7) An establishment shall inform MDA on the new lot number of the same batch, in order to issue new evaluation letter.

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Date: 19/11/2020

Order No.: OG 0117-2020

THIS IS TO CERTIFY THAT, ACCORDING TO THE COUNCIL DIRECTIVE 98/79/EC, OBELIS S.A. (O.E.A.R.C.) PERFORMED ALL NOTIFICATION DUTIES AND RESPONSIBILITIES AS THE EUROPEAN AUTHORIZED REPRESENTATIVE (EC REP) OF:

NAME: SHENZHEN LVSHIYUAN BIOTECHNOLOGY CO., LTD.

ADDRESS: 101, 201, 301, D BUILDING, NO. 2 INDUSTRIAL AVENUE, BUXIN VILLAGE, BUXIN COMMUNITY, DAPENG SUBDISTRICT OFFICE, DAPENG NEW DISTRICT, SHENZHEN, 518120, CHINA

AS STIPULATED AND DEMANDED BY THE AFOREMENTIONED DIRECTIVE.

The Manufacturer declares that the IVD devices comply with the Directive including all essential requirements.

The Manufacturer has provided Obelis s.a. (O.E.A.R.C.) with all the appropriate declarations according to the 98/79/EC Directive – article 10 requirements including the EC Declaration of Conformity confirming that his In-Vitro Diagnostics medical devices, as stipulated here above, are fulfilling the applicable requirements of the European Council Directive 98/79/EC

The notification of the following In-Vitro Diagnostic medical devices has been completed by Obelis s.a. (O.E.A.R.C.) on the 18/11/2020 in compliance with the European Council Directive 98/79/EC - article 10 requirements.

IN-VITRO DIAGNOSTIC MEDICAL DEVICES: PLEASE SEE ANNEX A - LIST OF DEVICES (1 PAGE, 1 DEVICE)

As of the 19/11/2020, and as long as the manufacturer will continue complying with the hereabove mentioned requirements* he therefore:

- Is required to affix the CE marking on these devices;
- Place these devices in the Territory of Belgium and/or the other EEA Member States (excluding territories not in alignment with Decision 2010/227/EU).

Obelis s.a. - O.E.A.R.C.
Registered Address:
Rue Cardinal Wauters 51
1030 Brussels

Mr. G. Elkayam CEO
Obelis sa

Obelis European Authorized Representative Center is a member of the European Association of Authorized Representatives (E.A.A.R.), ISO 9001 : 2015 and ISO 13485 : 2016 certified in accordance to the profession of a European Authorized Representative.

* This is not a CE mark and is only provided as a template for informational purposes.

**** This Certificate will be automatically void if the notification is rejected by the EU Authorities or upon termination of the EAR agreement.**


Order No.: OG 0117-2020

Ref No.: BS 0171-2020

Annex A - List of Devices

(Recital 29 of the Directive 98/79/EC on In Vitro Diagnostic Medical Devices)

#	Catalogue reference number	Commercial Name	Generic Device Term	Short description and intended use	GMDN/EDMS Code	Class
1	GF102B1	SARS-CoV-2 Antigen Rapid Test Kit	SARS-CoV-2 (COVID-19) Antigen Rapid Test Kit (Colloidal Gold)	SARS-CoV-2 (Covid-19) Antigen Rapid Test Kit is used to quantitatively detect antigen of the novel coronavirus in nasal swab.	15.04.80.19	Others

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Identification and Classification Report for Air Transport of Goods

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Beijing DGM Air Transport Technology Development Co., Ltd.

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依据联合国 GHS 制度第八修订版编/According to UN GHS (the 8th revised edition)

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ISSUED BY: ZHONGKE SERVICES OF TESTING (SHENZHEN) CO., LTD.

2022 年 01 月 印制 Jan. 2022 Print

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Applicant: Shenzhen Lvshiyuan Biotechnology Co., Ltd.